

Mactrac Pro Bed Screen Track System

Technical Data Sheet (TDS)

Product Description

The Mactrac Pro Bed Screen Track System is a premium quality and well-proven application for hospital bed screens, doctors' surgeries, rest homes, and other similar areas. The Mactrac Pro System uses extruded aluminium SC2 tracks with a structural rectangular hollow section with lipped glide channel below for a smooth glide. All tracks are available in powder coated Arctic White, Ivory, or Anodised Satin Silver track with Electro Natural Silver powder coated bracket. The thickness of the aluminium is a minimum of 2.1mm on the sidewalls and fixing walls to promote long life. The track can be fixed direct to the ceiling or suspended using the hanger tube system. Mactrac Pro can be bent to a 300mm radius and fabricated to a large variety of configurations. The versatile Mactrac Pro system can be fabricated to suit various layouts and situations.

Mactrac Pro Bed screen track SC2 Typical Layouts Typical Layouts

Typical Layout

Track Profile

Properties

Track Profile: 16.5mm x 25.4mm, available up to 6m lengths


Track Thickness: 2.1mm Standard Radius: 300mm

Glides: NRS1 polyethylene glide with stainless steel hook, NRS9 polyethylene glide with Acetyl Hook
Colour options: Arctic White, Ivory, Anodised Satin Silver with Electro Natural Silver powder coated brackets

Compatible with All SC2 fittings, anti-ligature, and quick-clip systems.

Componentry

Mactrac Pro is a complete system with a comprehensive range of componentry to ensure flexibility when installing.







Hospital Privacy Curtain Tracks

Brackets:	A
We have a range of Brackets; <i>wall</i> , <i>end</i> , and <i>swivel</i> brackets to ensure that the track can be installed, and work as expected in the most demanding of situations. Our Wall and End brackets retain the last glide and hook ensuring that the curtain stays close to the wall, maintaining patient privacy and dignity.	
Bed cubicle bays	2.
To facilitate the creation of bed bays, Mactrac provides a variety of track solutions. Mactrac offers both left and right-hand junction brackets, "T" brackets, and 4-way junction points.	
Suspended track, ceiling fixing:	
Tracks are often suspended from hospital ceilings for a variety of reasons including to allow air circulation and to reduce curtain fabric costs. The Mactrac hanger tube system provides a solution that can be tailored for a range of individual situations.	
Anti-ligature	
Mactrac offers two magnet-based systems. Both systems are triggered when weight in excess of over 20kgs is placed on the associated curtain. In the first system the track detaches from the bracket after 20kgs and in the "fuse" system, the "fuse" below the track above the curtain detaches. Rehang the curtain by simply putting the magnets together.	
Quick-clip	
The Quick-clip system is designed to accommodate a light disposable hospital curtain. It acts as an intermediary between the track and the curtain and facilitates coupling and decoupling from the track whilst using the well proven smooth and quiet Mactrac glides	

